


Photoelectric effect → phenomenon in which electrons are emitted from the surface of a metal when light strikes it.

⇒ There is a specific threshold freq., ν_0 , to emit an electron from a metal surface.



⇒ Below threshold freq., no e^- 's are emitted regardless of the intensity of the light.

⇒ For light $>$ freq., ν greater than threshold freq., the no. of e^- 's emitted \uparrow 's with the intensity of the light.

⇒ The KE of the emitted electrons \propto linearly \propto the freq. of the light used.

⇒ Threshold freq. is the minimum energy reqd to remove e^- from the metal surface $E_0 = h\nu_0$

If $\nu < \nu_0$ no e^- emitted
If $\nu > \nu_0$ e^- will emit

$$KE_{\text{electron}} = h\nu - h\nu_0$$

↓

KE of emitted electron

$$= E - E_0$$

Energy of incident photon Energy reqd to remove an electron from metal's surface

Photon electron spectroscopy :

is a technique in which a radⁿ $>$ enough energy can remove an e^- from an atom/molecule

when radⁿ is X-ray, then it is X-ray photoelectron spectroscopy

X-ray Photoelectron Spectroscopy

XPS ~~also~~ generates photoelectron spectrum of atoms/molecules.
~~Free Na-Atom~~

XPS also kld as ESCA (electron spectroscopy for chemical analysis) is the surface technique

→ provides both elemental & chemical state infⁿ of the material which can be analysed.

→ sample is ~~illustrated~~ ^{or irradiated} in vacuum by X-rays, and photoelectrons are emitted from the surface.

↓
K.E of emitted e⁻'s is characteristic of the element from which the photoelectron is originated.

The K.E of emitted photoelectron is equal to the diff. b/w the photon energy & binding energy of the electron.
 Φ = work function of spectrometer


$$KE = h\nu - BE - \Phi$$

This technique is inherently surface sensitive because the X-ray energy is low. The signal detected originates from the outer 1-10nm of a sample.

↓
K.E is measured by cylindrical mirror analyzer

↓
The spectrum plotted by the computer from the analyzer signal.

↓
B.E can be determined from the peak position & the elements present in the sample identified.



X-ray will emit
core^m electron

Components of XPS

- ① Source of X-rays \rightarrow emits photoelectron from the surface of sample
- ② An ultra high vacuum (UHV) \rightarrow to eliminate excessive surface contamination
- ③ An electron^{energy} analyser \rightarrow to measure k.E. of emitted e^- s
- ④ Magnetic field shielding
- ⑤ An electron detector system
- ⑥ A set of stage manipulators

Advantages \rightarrow measures composition of the surface
 \rightarrow measures chemical/electronic state of each element in the surface
 \rightarrow Uniformity of across the top surface

Applications : Polymer surface
in Industry : Corrosion

Adhesion


Thin film coating

Semiconductors

Disad - Expensive

Ultra vacuum req.

Slow



photoelectron
X-ray spectrophotometer

